


# CAPILLARY LINE

## APPLICATIONS

Irrigation of pots in a greenhouse or nursery



### STANDARD MODEL

Straight self-tapping stake.  
The standard stake accepts the  $\varnothing$  3.2 mm capillary tube



### COLIBRÌ MODEL

Elbow-shaped self-tapping stake.  
The colibrì stake accepts the  $\varnothing$  3.2 mm capillary tube  
Possibility of adjusting the angle of the capillary tube: 45° - 90°



### Codification

Code	Description
EU-RCP90885-1	Cut capillary tube 3.2 x 0.8 x 85 cm
EU-RCP908100-1	Cut capillary tube 3.2 x 0.8 x 100 cm
EU-RCP90800-X	Cut capillary tube 3.2 x 0.8 mm
EU-IPS0315	Colibrì stake (90°) for capillary tube $\varnothing$ 3.2 mm
EU-IPS0305-A	Stake for capillary tube from 2.9 to 3.2 mm, black
EU-IPS0305-A1	Stake for capillary tube from 2.9 to 3.2 mm, blue
EU-IPS0305-A2	Stake for capillary tube from 2.9 to 3.2 mm, red

